

Orbitrap Mass Analyser: a Tool for Titan Complex Molecular Content Exploration

Roland Thissen
Stefano Cornelli
Pascal Puget

Christelle Briois
Laurent Thirkell
Amirouche Boukrara
Jean-Pierre Lebreton

Patricia Beauchamp
Evan Neidholdt
Jack Beauchamp

Grenoble
Grenoble
Grenoble

Orléans
Orléans
Orléans
Orléans

Pasadena
Pasadena
Pasadena

Hervé Cottin
Noel Grand

Cyril Szopa
Nathalie Carrasco
Claude Chapuis
Abdel Bouabdellah

Cécile Engrand

Alexander Makarov




Creteil
Creteil
Versailles
Versailles
Versailles
Versailles

Orsay

Bremen
Toulouse
Paris


High mass-resolution for Titan?

Cassini (INMS) ionic densities in Titan ionosphere




	99.0109	resolution
C_7NH^+	99.0109	
C_8H_3^+	99.0235	7857
$\text{C}_4\text{N}_3\text{H}_9^+$	99.0797	1761
$\text{C}_5\text{N}_2\text{H}_{11}^+$	99.0923	7857
$\text{C}_6\text{NH}_{13}^+$	99.1049	7853
$\text{C}_7\text{H}_{15}^+$	99.1174	7920

Big molecule synthesis ...



Positive ions < 350



Negative ions > 1000


Bio molecules synthesis...

**Formation of amino acids and nucleotide bases
in a Titan atmosphere simulation experiment**

			P2CO	P2COi	P5CO	P5COi
	%N ₂	96.2	96.2	93.2	93.2	
	%CH ₄	2	2	5	5	
	%CO	1.8	1.8 C ¹⁸ O	1.8	1.8 C ¹⁸ O	
Name	Mass	Formula	Fig			
Nucleotide Base						
cytosine	111	C ₄ H ₅ N ₃ O	2	OT	OT	OT/GC-MS
uracil	112	C ₄ H ₄ N ₂ O ₂		OT	OT	OT/GC-MS
thymine	126	C ₅ H ₆ N ₂ O ₂		OT	OT	OT/GC-MS
adenine**	135	C ₅ H ₅ N ₅	2	OT		OT/GC-MS
guanine	151	C ₅ H ₅ N ₅ O		OT		OT/GC-MS
Biological Amino Acid						
glycine	75	C ₂ H ₅ NO ₂		OT		OT/GC-MS
alanine	89	C ₃ H ₇ NO ₂		OT		OT/GC-MS
serine	105	C ₃ H ₇ NO ₃		OT		OT
proline	115	C ₅ H ₉ NO ₂		OT		OT
valine	117	C ₅ H ₁₁ NO ₂		OT		OT
threonine	119	C ₄ H ₉ NO ₃		OT		
isoleucine/	131	C ₆ H ₁₃ NO ₂		OT		OT
leucine						
asparagine	132	C ₄ H ₈ N ₂ O ₃		OT		OT
glutamine	146	C ₅ H ₁₀ N ₂ O ₃		OT		OT
lysine	146	C ₆ H ₁₄ N ₂ O ₂		OT		OT
histidine	155	C ₆ H ₉ N ₃ O ₂	2	OT		OT
phenylalanine	165	C ₉ H ₁₁ NO ₂		OT		OT
arginine	174	C ₆ H ₁₄ N ₄ O ₂		OT		OT


Horst et al. Accepted in Astrobiology

Tholins mass spectrum at high resolution




Zoom on mass 262


070511_Lot1MeOH_b #663-697 RT: 19,32-20,28 Av. 35 NL: 1,02E5
T: FTMS + p ESI Full ms [50,00-500,00]



Isotopical abundances, example on complex mixture...



9 /




Orbitrap mass Analyser ?




Preliminary remark, what is orbitrap...

Mass Spectrometry :





Perfect quadro logarithmic potential



$$\omega = \sqrt{\frac{k}{m/z}}$$

record movement of
unperturbed ions during
 $>500\text{ms}$
pressure $<10^{-9} \text{ mbar!}$





How to inject and stabilize ions in this perfect potential? “Electrodynamic squeezing”

- a bunch of ions enters in the field, tangentially
- ions are squeezed towards central electrode by variation of its potential (deepening of well)
- The excitation of quadratic movement is inherent to “off axis” injection



Need to pulse ions in the trap!



Need for very stable HV potential
During measurement!

Orbitrap, potential

- # Ultra high resolution: 100 000 at mass 400, adjustable during mission, as it depends only on the integration time
- # Very small volume, lightweight : l=4 , f= 4 cm
- # Good detection Dynamic : 50 000 per spectrum
- # Positive or negative Ions as only one potential to invert
- # All ions are analysed simultaneously
- # No detector, no saturation, ...
- # no RF, no moving part
- # Ideal for solids or aerosols


Orbitrap for Titan




Sample handling

ionization



Ion accumulation




Study of optimized
MS instrument for the analysis of
aerosols/solids/liquids present in
the deep atmosphere/surface of
Titan



Orbitrap prototype in Orléans




Performances of prototype...




1 laser shot on Stainless steel

R=280 000



Tholins survive the laser desorption/ionization



Conclusions

- # Going back to Titan, we need to do much better than Cassini.
- # We know that very large (bio)organic molecules are produced/present in
 - Atmosphere **Aerosols**
 - Surface **Liquids / solids**
- # Request for very high resolution mass spectrometry to be characterized.
- # Orbitrap is a very promising analyser/detector for this purpose.
- # However big challenges remain:
 - Sample handling ?
 - Ionization method ?
 - High-Vacuum generation/maintenance ?
 - Beyond formulas, what about structural analysis ?
 - ...
- # Further info ? Come & visit posters **7A/B**
Briois, Cornelli and Thirkell